

			No	Type				No	Type				No	Type
<b>ADIN</b>						<b>BlueSkyBio</b>						<b>DIO</b>		
Touareg-S	55034	26	BIO Max NP Platform, All	55046	38	SM Submerged 3.8	55049	41						
Touareg Closefit UNP 2.75	55059	51	BIO Max DP Platform 3.0-3.5	55046	38	SM Submerged 5.0	55020	10						
			BIO Max DP Platform 4.3-5.0	55037	29									
<b>AIRES</b>						<b>Easy Implant</b>								
3.3	55049	41	BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	Ocean	55046	38						
			BIO Int Hex 3.5 Platform 4.7	55051	43	Lilas	55037	29						
			BIO Int Hex 4.5 Platform 4.7-7.0	55051	43									
<b>Alpha-Bio</b>														
Regular Platform	55034	26	BIO One-stage Narrow, Regular, Wide	55014	4									
Conical Platform	55065	57	BIO One-stage Short Regular & Wide	55014	4	<b>Edierre</b>								
			BIO Quattro Narrow 3.3	55061	53	3-lobe	55021	11						
			BIO Quattro Regular 4.1 - 5.6	55062	54	4-lobe	55025	16						
<b>Anker</b>														
SB III 3.5	55023	14	BIO Conus 12, 3.0 Platform 3.0	55064	56									
SB III 4.0-5.0	55020	10	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38	<b>Global D</b>								
			BIO Conus 12, 4.5 x 6mm	55046	38	In-Kone 4.0-5.0	55024	15						
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7									
<b>Ankylos</b>						<b>gmi</b>								
C/X	55025	16	BIO Conus 12, 5.0 x 6mm	55046	38	avantgard RP 3.75-4.25	55066	58						
			<b>Bredent</b>			avantgard WP 4.75-5.75	55075	67						
<b>Anthogyr</b>			SKY			55025	16	frontier RP,WP	55051	43				
Axiom, all	55065	57	copaSKY	55027	19	phoenix NP 3.3	55038	30						
			mini2SKY	55033	25	phoenix Std 3.3 - Wide 5.0	55011	1						
<b>Argon</b>			fast & fixed Abutment			55083	75							
K3pro Rapid 3.0 - 3.5	55026	17				<b>Hahn</b>								
K3pro Rapid 4.0 - 6.0	55031	23	<b>BTI</b>			Tapered 3.0			55058	50				
K3pro Rapid Platf 3mm Short	55044	36	Interna Universal	55051	43	Tapered 3.5	55046	38						
			Interna Wide	55051	43	Tapered 4.3	55046	38						
<b>Astra Tech</b>						Tapered 5.0	55037	29						
Osseospeed EV 3.0	55058	50	<b>BTK</b>			Tapered 7.0	55048	40						
Osseospeed EV 3.6	55046	38	External Hex EN	55038	30									
Osseospeed EV 4.2	55057	49	External Hex ER	55011	1	<b>Hiossen</b>								
Osseospeed EV 4.8	55055	47	External Hex EW	55011	1	HS System	55028	20						
Osseospeed EV 5.4	55048	40	Internal Hex IR	55049	41	ET System	55037	29						
Osseospeed TX 3.0 S	55064	56	Internal Hex IM	55051	43	ET 3.2	55047	39						
Osseospeed TX 3.5-4.0 S	55046	38	Internal Hex IW	55051	43	ET mini	55015	5						
Osseospeed TX 4.5-5.0 S	55017	7	Hexagonal Morse Taper KR	55081	73									
			Hexagonal Morse Taper KW	55056	48	<b>IDI Evolution</b>								
			Octagonal Morse Taper SR	55014	4	Standard implant 3.75-4.0	55004	18						
<b>Avinent</b>														
External, Coral & Ocean 3.5	55038	30				<b>Implance</b>								
External, Coral & Ocean 4.1	55011	1	<b>Camlog</b>			Tissue level			55014	4				
External, Coral 5.1	55013	3	Screw-Line, Root-Line 3.3	55036	28	Bone level Regular	55037	29						
Internal, Coral & Ocean 3.5, 4.1	55051	43	Screw-Line, Root-Line 3.8, 4.3	55031	23	Bone level 3.3	55065	57						
Conical, Ocean 3.5	55066	58	Screw-Line, Root-Line 5.0, 6.0	55032	24									
Conical, Ocean 4.1	55004	18	iSy	55065	57	<b>Implant Direct</b>								
<b>BEGO</b>			<b>Conelog</b>			Legacy 2&3&4 3.0/3.2			55012	2				
Semados Mini	55012	2	Screw-Line 3.3	55066	58	Legacy 2&3&4 3.5/3.7 - 5.7/7.0	55051	43						
Semados S / RI 3.25-5.5	55034	26	Screw-Line 3.8, 4.3	55066	58									
Semados RS / RSX 3.0	55082	74	Screw-Line 5.0	55075	67	<b>Implant Swiss</b>								
Semados RS / RSX 3.75-5.5	55034	26				Tissue level	55014	4						
Semados SC / SCX 3.25-5.5	55034	26	<b>Dentack</b>			Bone level Regular	55037	29						
			Internal Hex	55056	48	Bone level 3.3	55065	57						
<b>BioHorizons</b>														
Laser-Lok 3.0	55034	26	<b>Dentegris</b>			<b>Implantvision</b>								
Internal Laser-Lok 3.5 Platform	55049	41	3.3	55034	26	Internal Hex 3.7-5.7	55034	26						
Internal Laser-Lok 4.5 Platform	55051	43	3.75	55057	49									
Internal Laser-Lok 5.7 Platform	55051	43	4.5-5.5	55049	41	<b>Intra-lock</b>								
External 4.0 Platform	55011	1				Intrahex 3.4 - 4.5	55051	43						
			<b>C-Tech</b>			Gold & Blue 3.4 - 6.0	55051	43						
<b>Biomet/3i</b>			EL			55067	59	InDex 3.4 - 3.75	55015	5				
Certain internal Purple	55036	28				InDex 4.0 - 4.3	55020	10						
Certain internal Purple Platform 3.4	55066	58	<b>Dentium</b>											
Certain internal Blue, Yellow, Green	55066	58	Superline	55037	29	<b>J Dental Care</b>								
			NR Line	55066	58	JDEvolution 3.7-6.0	55067	59						
<b>Bionnovation</b>			Simple Line II GH 1,6			55014	4							
Biomorse	55023	14				<b>Keystone</b>								
<b>Biotech Dental</b>			<b>DEVA</b>			Prima / Genesis Narrow			55049	41				
Kontakt, Kontakt-S 3.6 - 4.8	55045	37	Bone Level 3.5	55064	56	Prima / Genesis Regular	55067	59						
			Bone Level 4.5	55051	43	Prima / Genesis Wide	55041	33						
						Restore Narrow	55079	71						
						Restore Regular	55011	1						
						Restore Wide	55013	3						

	No	Type		No	Type		No	Type
<b>Klockner</b>			<b>OCO Biomedical</b>			<b>Sweden &amp; Martina</b>		
VEGA Micro 3.0	55074	66	Engage 3.25	55027	19	Syra, Syra SP 3.9-5.0	55011	1
VEGA Narrow 3.5	55065	57	Engage 4.0-6.0	55051	43	Outlink <sup>2</sup> 3.75, 4.10, 5.0	55011	1
VEGA Regular 4.0, 4.5	55034	26	ERI 3.25	55069	61	Shelta 3.8-5.0	55025	16
KL Narrow 3.3	55012	2	ERI 4.0-5.0	55014	4	Premium Straight & SP 3.3-3.8	55025	16
KL Regular 3.7, 4.2	55011	1	TSI 3.25	55068	60	Kohno Straight & SP 3.3-3.8	55025	16
KL Wide 4.7	55013	3	TSI 4.0-5.0	55011	1	Premium Straight & SP 4.25-6.0	55004	18
SK2	55038	30				Kohno Straight & SP 4.25-6.0	55004	18
Essential Cone	55067	59	<b>Osstem</b>			Prama 3.8 - 5.0	55079	71
S3M 3.1	55079	71	US Regular	55011	1			
			US Mini	55012	2	<b>Tag Dental</b>		
<b>Konmet</b>			US Wide	55013	3	Axis	55051	43
Classic 4.0, 4.8	55014	4	US Wide PS	55011	1	massif	55051	43
Cylindrical 3.5-4.0	55065	57	SS Regular, Wide	55028	20			
Cylindrical 4.5-5.5	55057	49	TS Regular, Reg Short, Diam 6 & 7	55037	29	<b>Thommen</b>		
Short Cylindrical 5.0-6.0	55078	70	TS Mini 3.5	55015	5	SPI 3.5	55021	11
Conical 4.4/3.4	55065	57	TS Mini 3.0	55047	39	SPI 4.0	55065	57
Conical 4.8/3.8, 5.5/4.5	55057	49				SPI 4.5-6.0	55066	58
			<b>Paltop</b>					
<b>Lasak</b>			Standard 3.75-5.0	55034	26	<b>TRI Dental implants</b>		
BioniQ QR	55066	58	Narrow	55027	19	Tissue level Octa	55014	4
BioniQ QN	55021	11	PCA Conical	55066	58	Tri-Vent, all	55034	26
						Narrow	55025	16
<b>Leader</b>			<b>Phibo</b>					
Internal Hex, 3.3 - 5.5 mm	55034	26	Aurea EVO NP	55027	19	<b>Xive</b>		
Cono-In	55067	59	Aurea EVO RP, WP	55051	43	S-Plus 3.0	55030	22
			TSA S3	55074	66	S-Plus 3.4	55061	53
<b>Maxillent</b>			TSA S4, S5	55034	26	S-Plus 3.8-5.5	55062	54
iRaise	55049	41	TSH S3, S4, S5	55011	1	TG Plus 3.4-4.5	55018	8
<b>Medentika</b>			<b>Prodent italia</b>			<b>Zimmer</b>		
Microcone NI	55021	11	Prime, Twinner: Orange Platform 3.3	55027	19	Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Microcone RI	55081	73	Prime, Twinner: Fuchsia Platform 3.6	55051	43	Screw-Vent, Tapered Screw-Vent 4.7	55051	43
Quattrocone RI	55081	73	Prime, Twinner: Green Platform 4.0	55051	43			
			Prime, Twinner: Yellow Platform 4.5	55020	10	<b>Zuga</b>		
<b>Megagen</b>			Prime: Blue Platform 5.0	55020	10	All	55027	19
Anyridge	55051	43						
			<b>Resista</b>					
<b>MIS</b>			WIS 3.0-5.0	55049	41			
Conical connection Narrow (C1, V3)	55081	73						
Conical connection Standard (C1, V3)	55067	59	<b>Schutz Dental</b>					
Conical connection Wide (C1)	55041	33	Impla cylindrical	55052	44			
Seven and M4 Narrow	55027	19						
Seven and M4 Standard	55049	41	<b>SGS Implant</b>					
Seven Wide 5.0	55051	43	P1, P7 Platform	55049	41			
			P1D, P7D	55049	41			
<b>Neo Biotech</b>			P7N	55059	51			
IS II 4.0	55020	10						
IS III 4.0	55020	10	<b>SIC</b>					
			Ace 3.4 - 5.0	55052	44			
<b>Neodent</b>			Vantage 4.2-5.2	55054	46			
Drive CM	55025	16						
Titamax WS Cortical	55056	48	<b>Southern Implants</b>					
			Tri-Nex 4.3-5.0	55010	13			
<b>Neoss</b>			Tri-Nex 3.5	55022	12			
Diameter 3.5-5.5	55004	18	External Hex	55011	1			
Diameter 3.25	55027	19	MAX Ext Hex	55011	1			
Access Abutment	55060	52	Co-Axis Ext Hex	55011	1			
			Piccolo Ext Hex	55050	42			
<b>Nobel Biocare</b>			IT-Connection	55014	4			
Brånemark System Ext hex RP	55011	1	Deep Conical 3.0	55064	56			
Brånemark System Ext hex NP	55012	2	Deep Conical 3.5-4.0	55046	38			
Brånemark System Ext hex WP	55013	3						
Replace Select Tri-lobe NP	55022	12	<b>Spiral Tech</b>					
Replace Select Tri-lobe RP, WP	55010	13	ESi, all	55051	43			
Replace Select Tri-lobe 6.0	55032	24	Ultimate	55051	43			
Nobel Active 3.0	55058	50						
Nobel Active Narrow	55046	38	<b>Straumann</b>					
Nobel Active Regular	55037	29	Tissue level Synocta	55014	4			
Nobel Active Wide 5.5	55048	40	Tissue level NNC	55076	68			
MUA NP & RP	55080	72	Bone level / BLT NC	55061	53			
			Bone level / BLT RC	55062	54			
<b>NucleOss</b>			Bone Level /BLT SC	55064	56			
Tpure Bone Level	55074	66						
T6 Bonel Level	55004	18						